



Titan™ 80-300

Expanding the boundaries
and achieving new results





spherical aberration can be corrected and the typical C_s constant is of less importance: the point resolution now equals the information limit. The information limit is governed by the stability of the microscope and is the best measure for corrected microscopy quality and performance. The requirement for unsurpassed stability calls for a corrected microscope design that complies with the stringent needs for maximum mechanical, electronic and thermal stability, as well as precision alignment of the advanced components.

The Titan 80-300™ microscope incorporates a newly designed platform dedicated to the principles of ultimate stability, ultimate performance and ultimate flexibility for corrector and monochromator technology and its applications. The microscope transfers information deep into sub-Ångström resolution making way for the highest performance available in both TEM and STEM. Titan obtains a better lateral resolution than 1 Ångström (see figure below) and energy resolution down to 0.1 eV. This presents new information of the electronic properties of materials such as bonding states or band gaps with unprecedented spatial resolution. An innovative modular and patented design makes a field upgrade with a C_s corrector possible, allowing you a two step approach to corrected microscopy in your laboratory.


Conquering the sub-Ångström Era

Nanoresearch strives to increase our understanding of nanostructures and functional materials by linking the macroscopic material properties to the properties at the smallest level of detail: down to atoms. This continuing drive pushes for the better understanding of local atom organization and for the characterization of atoms, chemical bonding and even the electronic structure.


The recognized and ever growing need to truly image and analyze nanostructures with the highest level of detail implies setting new standards for achieving new results in nanoresearch microscopy. Sub-Ångström aberration-corrected analytical microscopy is the benchmark and presents a new set of rules in this new era. A classical electron microscope's fundamental

The stability, performance and ease-of-use of Titan enables corrected microscopy to be taken to the next level where new discoveries on the structure-property relationships of materials become possible at ever-decreasing scales. Titan is poised to bring electron microscopy into the new era by expanding the boundaries and achieving new results in nanoresearch.

Titan with monochromator




Energy resolution at 300kV using FEI proprietary monochromator with embedded Gatan Tridiem 866 is better than 0.2 eV with 1 second exposure time.




High-resolution EELS spectrum of GaAs and GaAlAs using Titan's monochromator technology. Uniquely, after subtraction of the zero-loss peak, bandgap thresholds can be measured as well as the energy shift due to the Al doping.

Ultimate performance, ultimate stability, ultimate flexibility

Titan basic performance (no correctors)




HR STEM. Grain boundary in gold in $\langle 110 \rangle$ direction. (raw data). Missing atomic columns are visible at the boundary
Sample courtesy C. Kisielowski NCEM Berkley (USA).




HR TEM. New Rules: Classical " C_s " is now irrelevant and point resolution now equals information limit. Information limit is all about stability. Titan's Information Limit: Young's Fringes of polycrystalline gold. Fringes extend below 1.0 \AA .


Titan with C_s corrector



Atomic resolution HR-STEM image of diamond $\langle 110 \rangle$ using Titan's C_s probe corrector. Diamond dumbbells are shown at a spacing of 0.9 \AA . The twin and grain boundaries can be clearly imaged.




Atomic resolution HR-TEM direct image of Pt/Co nano particles using Titan's C_s image corrector. Individual atom columns and particle boundaries can be clearly resolved.




Atomic resolution HR-STEM image of Ge $\langle 112 \rangle$ using Titan's C_s probe corrector. The Germanium atomic distance of 0.81 \AA can be resolved in the real image (left). In Fourier space resolutions down to 0.63 \AA are resolved (right).

Courtesy: Ernst Ruska Center, Germany.



Atomic resolution HR-TEM direct image of Silicon $\langle 110 \rangle$ using Titan's C_s image corrector. Individual atom columns can be clearly resolved.

Elemental maps using Titan's simultaneous acquisitions of EELS and EDS signals in STEM. Elemental maps reveal the entire chemistry (low Z and high Z) of the device under investigation. 200×200 pixels map with a scanning probe smaller than 5 \AA .



FEI Company

World Headquarters
and North American Sales
5350 NE Dawson Creek Drive
Hillsboro, OR 97124-5793 USA
Tel: +1 503 726 7500
Fax: +1 503 726 7509

European Sales
Tel: +31 40 23 56 110
Fax: +31 40 23 56 612

Asia-Pacific Sales
Tel: +65 6272 0050
Fax: +65 6272 0034

Japan Sales
Tel: +81 3 3740 0970
Fax: +81 3 3740 0975

e-mail: sales@fei.com
www.fei.com

©2006. We are constantly improving the performance of our products, so all specifications are subject to change without notice. The FEI logo, Tools for Nanotech and Titan are trademarks of FEI Company. Windows is a trademark of Microsoft Corporation. Tridien is a product of Gatan Incorporated.

TÜV Certification for design, manufacture, installation and support of focused ion- and electron-beam microscopes for the NanoElectronics, NanoResearch and NanoBiology markets.

